

Pipet Tip Compatibility

Pipette Brand		Axygen Axyjet (Single channel)					Axygen Axyjet (Multichannel)				Nichiyu/Oxford Benchmate			Biohit Proline (Single channel)					
Model / Range		0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-50µl	0.2-10µl	5-100µl	10-250µl	10-500µl	50-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓					✓				✓			✓			✓	
	21-0010	✓	✓					✓				✓			✓			✓	
1-200ul	20-0200			✓	✓	✓	✓		✓	✓	✓			✓		✓	✓		
	21-0200			✓	✓	✓	✓		✓	✓	✓			✓		✓	✓		
100-1000ul	20-1000																	✓	✓
	21-1000																	✓	✓
100-1250ul	20-1250																	✓	✓
	20-1251																	✓	✓
	21-1250																	✓	✓

Pipette Brand		Biohit Proline (Multichannel)			Biohit e-line (Single channel)				Biohit e-line (Multi channel)			Biohit m-line (Single channel)			Biohit m-line (Multichannel)				
Model / Range		0.2-10µl	5-100µl	25-250µl	1-10µl	5-120µl	10-300µl	50-1000µl	0.2-10µl	5-120µl	10-300µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	5-100µl	30-300µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓			✓				✓									✓	
	21-0010	✓			✓				✓									✓	
1-200ul	20-0200		✓	✓		✓	✓			✓	✓		✓	✓	✓			✓	✓
	21-0200		✓	✓		✓	✓			✓	✓		✓	✓	✓			✓	✓
100-1000ul	20-1000							✓								✓			
	21-1000							✓								✓			
100-1250ul	20-1250							✓								✓			
	20-1251							✓								✓			
	21-1250							✓								✓			

Pipette Brand		Capp Aero (Multichannel)							BrandTech Transferpette (Single channel)							MatrixImpact/Impact(s/mc)			
Model / Range		0.2-2µl	0.5-10µl	5-50µl	10-100µl	25-200µl	50-300µl	100-1000µl	0.1-1µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	25-250µl	100-1000µl	20-200µl	5-100µl	30-300µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓						✓	✓									
	21-0010	✓	✓						✓	✓									
1-200ul	20-0200			✓	✓	✓	✓				✓	✓	✓	✓		✓			
	21-0200			✓	✓	✓	✓				✓	✓	✓	✓		✓			
100-1000ul	20-1000							✓							✓				
	21-1000							✓							✓				
100-1250ul	20-1250							✓							✓				✓
	20-1251							✓							✓				✓
	21-1250							✓							✓				✓

Pipet Tip Compatibility

Pipette Brand		BrandTech Transferpette (Multichaanel)						Rainin EDP3 (Single channel)						Rainin Latch Mode					
Model / Range		0.5-10µl	2-20µl	2.5-25µl	5-50µl	10-100µl	20-200µl	30-300µl	E-MIC-10, 0.5-10µl	E2-25, 2.5-25µl	E2-100, 10-100µl	E2-250, 25-250µl	E2-1000, 100-1000µl	R2, 0.1-2µl	R10, 0.2-10µl	R20, 2-20µl	R100, 10-100µl	R200, 20-200µl	R1000, 100-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓						✓					✓	✓				
	21-0010	✓	✓						✓					✓	✓				
1-200ul	20-0200			✓	✓	✓	✓		✓	✓	✓				✓		✓	✓	
	21-0200			✓	✓	✓	✓		✓	✓	✓				✓		✓	✓	
100-1000ul	20-1000							✓					✓						✓
	21-1000							✓					✓						✓
100-1250ul	20-1250							✓					✓						✓
	20-1251							✓					✓						✓
	21-1250							✓					✓						✓

Pipette Brand		Labnet Labpette (Single chaanel)					Labnet Labpette (Multichaanel)				Labnet Biopette E (Single chaanel)			Labnet Biopette E (Multichaanel)					
Model / Range		0.1-2µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-200µl	100-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓					✓				✓	✓			✓	✓		
	21-0010	✓	✓					✓				✓	✓			✓	✓		
1-200ul	20-0200			✓	✓	✓			✓	✓	✓		✓					✓	
	21-0200			✓	✓	✓			✓	✓	✓		✓					✓	
100-1000ul	20-1000						✓						✓						✓
	21-1000						✓						✓						✓
100-1250ul	20-1250						✓						✓						✓
	20-1251						✓						✓						✓
	21-1250						✓						✓						✓

Pipette Brand		Eppendorf Electronic				Socorex Calibra/Accura (Single chaanel)						VWR AND VWR Signature(Multichannel)				Rainin EDP3 (Multichannel)				
Model / Range		0.5-10µl	5-100µl	2-200µl	50-1000µl	0.1-2µl	0.2-2µl	0.5-10µl	1-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	EP-M8-10, 0.5-10µl	EP-M8250, 10-250µl	P200-M8, 20-200µl
Volume	Product No.	Universal Fit Filter Tips																		
0.5-10ul	20-0010	✓				✓	✓	✓	✓					✓				✓		
	21-0010	✓				✓	✓	✓	✓					✓				✓		
1-200ul	20-0200		✓	✓						✓	✓	✓			✓	✓	✓		✓	✓
	21-0200		✓	✓						✓	✓	✓			✓	✓	✓		✓	✓
100-1000ul	20-1000				✓								✓							✓
	21-1000				✓								✓							✓
100-1250ul	20-1250				✓								✓							✓
	20-1251				✓								✓							✓
	21-1250				✓								✓							✓

Pipet Tip Compatibility

Pipette Brand		VWR AND VWR Signature (Single channel)								Socorex Calibra/Accura (Multichannel)					EppendorfReference (single channel) (multichannel)				
Model / Range		0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-250µl	100-1000µl	1-10µl	5-50µl	10-100µl	20-200µl	50-200µl	0.1-10µl	2-20µl	10-100µl	20-200µl	100-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓							✓					✓				
	21-0010	✓	✓							✓					✓				
1-200ul	20-0200			✓	✓	✓	✓	✓			✓	✓	✓	✓		✓	✓	✓	
	21-0200			✓	✓	✓	✓	✓			✓	✓	✓	✓		✓	✓	✓	
100-1000ul	20-1000								✓						✓	✓	✓	✓	
	21-1000								✓						✓	✓	✓	✓	
100-1250ul	20-1250								✓										✓
	20-1251								✓										✓
	21-1250								✓										✓

Pipette Brand		Eppendorf Research (singlechannel)						Eppendorf Research (multichannel)						Thermo Scientif Finnpiquette (single channel)					
Model / Range		0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	25-200µl	100-1000µl	0.5-10µl	5-100µl	50-100µl	20-300µl	30-300µl	50-1200µl	0.2-10µl	0.5-10µl	5-40µl	40-200µl	100-1000µl	200-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓					✓						✓	✓				
	21-0010	✓	✓					✓						✓	✓				
1-200ul	20-0200				✓	✓				✓	✓	✓	✓		✓	✓	✓		
	21-0200				✓	✓				✓	✓	✓	✓		✓	✓	✓		
100-1000ul	20-1000				✓	✓				✓	✓	✓	✓		✓	✓		✓	✓
	21-1000				✓	✓				✓	✓	✓	✓		✓	✓		✓	✓
100-1250ul	20-1250						✓						✓					✓	✓
	20-1251						✓						✓					✓	✓
	21-1250						✓						✓					✓	✓

Pipette Brand		Thermo Scientif Finnpiquette (multichannel)		Gilson Pipetma (single channel)					
Model / Range		5-50µl	50-300µl	P2,0.1-2µl	P10, 0.2-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl
Volume	Product No.	Universal Fit Filter Tips							
0.5-10ul	20-0010			✓	✓				
	21-0010			✓	✓				
1-200ul	20-0200	✓	✓			✓	✓	✓	
	21-0200	✓	✓			✓	✓	✓	
100-1000ul	20-1000	✓	✓			✓	✓	✓	
	21-1000								✓
100-1250ul	20-1250								✓
	20-1251								✓
	21-1250								✓

Pipet Tip Compatibility

Pipette Brand		Eppendorf Research Plus						Gilson Pipetman(PG)					Gilson Pipetman(PN)					Gilson Pipetman(PL)								
Model / Range		0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100,10-100µl	P200, 20-200µl	P1000, 100-1000µl	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl	
Volume	Product No.	Universal Fit Filter Tips																								
0.5-10ul	22-0010S	✓	✓					✓	✓					✓	✓					✓	✓					
	22-0011S	✓	✓					✓	✓					✓	✓					✓	✓					
1-20ul	22-0020S			✓					✓							✓							✓			
	22-0021S			✓					✓							✓							✓			
1-50ul	22-0050S				✓						✓	✓					✓	✓					✓	✓		
	22-0051S				✓						✓	✓					✓	✓					✓	✓		
1-100ul	22-0100S				✓						✓	✓					✓	✓					✓	✓		
	22-0101S				✓						✓	✓					✓	✓					✓	✓		
1-300ul	22-0200S					✓						✓						✓						✓		
	22-0201S					✓						✓						✓						✓		
100-1000ul	22-1000S						✓						✓						✓					✓		
	22-1001S						✓						✓						✓					✓		

Pipette Brand		Thermo Finnpipette (F1)										Thermo Finnpipette (F2)										Thermo Finnpipette (F3)														
Model / Range		0.2-2µl	1-10µl	1-10µl(µTip)	2-20µl	2-20µl(µTip)	5-50µl	5-50µl(µTip)	10-100µl	20-200µl	30-300µl	100-1000µl	0.2-2µl(µTip)	1-10µl	1-10µl(µTip)	2-20µl	2-20µl(µTip)	5-50µl	5-50µl(µTip)	10-100µl	20-200µl	100-1000µl	0.2-2µl	1-10µl	1-10µl(µTip)	2-20µl	2-20µl(µTip)	5-50µl	5-50µl(µTip)	10-100µl	20-200µl	100-1000µl				
Volume	Product No.	Universal Fit Filter Tips																																		
0.5-10ul	22-0010S	✓	✓										✓										✓	✓												
	22-0011S	✓	✓										✓										✓	✓												
1-20ul	22-0020S			✓											✓											✓										
	22-0021S			✓											✓											✓										
1-50ul	22-0050S				✓			✓	✓								✓										✓	✓								
	22-0051S				✓			✓	✓								✓										✓	✓								
1-100ul	22-0100S				✓			✓	✓								✓										✓	✓								
	22-0101S				✓			✓	✓								✓										✓	✓								
1-300ul	22-0200S					✓				✓																		✓	✓							
	22-0201S					✓				✓																		✓	✓							
100-1000ul	22-1000S										✓																	✓								
	22-1001S										✓																	✓								

Pipet Tip Compatibility

Pipette Brand		Thermo Finnpiipette (F1)								Thermo Finnpiipette (F2)								Thermo Finnpiipette (F3)									
Model / Range		0.1-2.5µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-300µl	100-1000µl	0.1-3µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	30-300µl	100-1000µl	0.1-1µl	0.1-2.5µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	30-300µl	100-1000µl	
Volume	Product No.	Universal Fit Filter Tips																									
0.5-10ul	22-0010S	✓	✓							✓	✓							✓	✓	✓							
	22-0011S	✓	✓							✓	✓								✓	✓	✓						
1-20ul	22-0020S			✓								✓									✓						
	22-0021S			✓								✓									✓						
1-50ul	22-0050S				✓	✓	✓	✓					✓	✓	✓	✓						✓	✓	✓	✓	✓	✓
	22-0051S				✓	✓	✓	✓					✓	✓	✓	✓						✓	✓	✓	✓	✓	✓
1-100ul	22-0100S					✓	✓	✓					✓	✓	✓	✓							✓	✓	✓	✓	✓
	22-0101S					✓	✓	✓					✓	✓	✓	✓							✓	✓	✓	✓	✓
1-300ul	22-0200S						✓	✓						✓	✓	✓								✓	✓	✓	✓
	22-0201S						✓	✓						✓	✓	✓								✓	✓	✓	✓
100-1000ul	22-1000S								✓								✓										✓
	22-1001S								✓								✓										✓

Pipette Brand		Eppendorf Research Plus						Gilson Pipetman(PG)					Gilson Pipetman(PN)					Gilson Pipetman(PL)									
Model / Range		0.1-2.5ul	0.5-10ul	2-20ul	10-100ul	20-200ul	100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul		
Volume	Product No.	Universal Fit Filter Tips																									
0.5-10ul	23-0010S	✓	✓					✓	✓					✓	✓					✓	✓						
	23-0011S	✓	✓					✓	✓					✓	✓					✓	✓						
1-20ul	23-0020S			✓						✓						✓							✓	✓			
	23-0021S			✓						✓						✓							✓	✓			
1-50ul	23-0050S				✓	✓	✓				✓	✓				✓	✓	✓					✓	✓	✓	✓	✓
	23-0051S				✓	✓	✓				✓	✓				✓	✓	✓					✓	✓	✓	✓	✓
1-100ul	23-0100S				✓	✓	✓				✓	✓				✓	✓	✓					✓	✓	✓	✓	✓
	23-0101S				✓	✓	✓				✓	✓				✓	✓	✓					✓	✓	✓	✓	✓
1-300ul	23-0200S					✓	✓					✓						✓	✓					✓	✓	✓	✓
	23-0201S					✓	✓					✓						✓	✓					✓	✓	✓	✓
100-1000ul	23-1000S						✓					✓						✓	✓							✓	✓
	23-1001S						✓					✓						✓	✓							✓	✓

Pipet Tip Compatibility

Pipette Brand		Thermo Finnpiquette (F1)										Thermo Finnpiquette (F2)								Thermo Finnpiquette (F3)															
Model / Range		0.2-2ul	1-10ul	1-10ul(μTip)	2-20ul	2-20ul(μTip)	5-50ul	5-50ul(μTip)	10-100ul	20-200ul	30-300ul	100-1000ul	0.2-2ul(μTip)	1-10ul	1-10ul(μTip)	2-20ul	2-20ul(μTip)	5-50ul	5-50ul(μTip)	10-100ul	20-200ul	100-1000ul	0.2-2ul	1-10ul	1-10ul(μTip)	2-20ul	2-20ul(μTip)	5-50ul	5-50ul(μTip)	10-100ul	20-200ul	100-1000ul			
Volume	Product No.	Universal Fit Filter Tips																																	
0.5-10ul	23-0010S	✓	✓										✓										✓	✓											
	23-0011S	✓	✓										✓										✓	✓											
1-20ul	23-0020S				✓										✓											✓									
	23-0021S				✓										✓											✓									
1-50ul	23-0050S					✓			✓	✓							✓											✓							
	23-0051S					✓			✓	✓							✓											✓							
1-100ul	23-0100S					✓			✓	✓							✓											✓		✓					
	23-0101S					✓			✓	✓							✓											✓		✓					
1-300ul	23-0200S										✓																			✓					
	23-0201S										✓																			✓					
100-1000ul	23-1000S										✓												✓									✓			
	23-1001S										✓												✓									✓			

Pipette Brand		BIOHIT Proline								BIOHIT Proline (Plus)								BRAND																	
Model / Range		0.1-2.5ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	50-300ul	100-1000ul	0.1-3ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	30-300ul	100-1000ul	0.1-1ul	0.1-2.5ul	0.5-10ul	2-20ul	5-50ul	10-100ul	20-200ul	30-300ul	100-1000ul									
Volume	Product No.	Universal Fit Filter Tips																																	
0.5-10ul	23-0010S	✓	✓							✓	✓							✓	✓	✓															
	23-0011S	✓	✓							✓	✓							✓	✓	✓															
1-20ul	23-0020S			✓								✓										✓													
	23-0021S			✓								✓										✓													
1-50ul	23-0050S				✓	✓	✓	✓				✓	✓	✓	✓							✓		✓	✓	✓									
	23-0051S				✓	✓	✓	✓				✓	✓	✓	✓							✓		✓	✓	✓									
1-100ul	23-0100S					✓	✓	✓					✓	✓	✓									✓	✓	✓									
	23-0101S					✓	✓	✓					✓	✓	✓									✓	✓	✓									
1-300ul	23-0200S						✓	✓						✓	✓										✓	✓									
	23-0201S						✓	✓						✓	✓										✓	✓									
100-1000ul	23-1000S								✓																				✓			✓			
	23-1001S								✓																				✓			✓			