

# Pipet Tip Compatibility

| Pipette Brand |             | Axygen Axyet (Single channel) |          |        |        |          |          | Axygen Axyet (Multichannel) |        |          |          | Nichiyu/Oxford Benchmate |        |         | Biohit Proline (Single channel) |         |          |          |           |
|---------------|-------------|-------------------------------|----------|--------|--------|----------|----------|-----------------------------|--------|----------|----------|--------------------------|--------|---------|---------------------------------|---------|----------|----------|-----------|
| Model / Range |             | 0.1-2µl                       | 0.5-10µl | 2-20µl | 5-50µl | 10-100µl | 20-200µl | 1-10µl                      | 5-50µl | 20-200µl | 50-300µl | 0.5-10µl                 | 2-20µl | 10-50µl | 0.2-10µl                        | 5-100µl | 10-250µl | 10-500µl | 50-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips     |          |        |        |          |          |                             |        |          |          |                          |        |         |                                 |         |          |          |           |
| 0.5-10ul      | 20-0010     | ✓                             | ✓        |        |        |          |          | ✓                           |        |          |          | ✓                        |        |         | ✓                               |         |          | ✓        |           |
|               | 21-0010     | ✓                             | ✓        |        |        |          |          | ✓                           |        |          |          | ✓                        |        |         | ✓                               |         |          | ✓        |           |
| 1-200ul       | 20-0200     |                               |          | ✓      | ✓      | ✓        | ✓        |                             | ✓      | ✓        | ✓        |                          | ✓      | ✓       |                                 | ✓       | ✓        |          |           |
|               | 21-0200     |                               |          | ✓      | ✓      | ✓        | ✓        |                             | ✓      | ✓        | ✓        |                          | ✓      | ✓       |                                 | ✓       | ✓        |          |           |
| 100-1000ul    | 20-1000     |                               |          |        |        |          |          |                             |        |          |          |                          |        |         |                                 |         |          | ✓        | ✓         |
|               | 21-1000     |                               |          |        |        |          |          |                             |        |          |          |                          |        |         |                                 |         |          | ✓        | ✓         |
| 100-1250ul    | 20-1250     |                               |          |        |        |          |          |                             |        |          |          |                          |        |         |                                 |         |          | ✓        | ✓         |
|               | 20-1251     |                               |          |        |        |          |          |                             |        |          |          |                          |        |         |                                 |         |          | ✓        | ✓         |
|               | 21-1250     |                               |          |        |        |          |          |                             |        |          |          |                          |        |         |                                 |         |          | ✓        | ✓         |

| Pipette Brand |             | Biohit Proline<br>(Multichannel) |         |          | Biohit e-line<br>(Single channel) |         |          |           | Biohit e-line<br>(Multi chaanel) |         |          | Biohit m-line<br>(Single channel) |        |          | Biohit m-line<br>(Multichannel) |            |          |         |          |
|---------------|-------------|----------------------------------|---------|----------|-----------------------------------|---------|----------|-----------|----------------------------------|---------|----------|-----------------------------------|--------|----------|---------------------------------|------------|----------|---------|----------|
| Model / Range |             | 0.2-10µl                         | 5-100µl | 25-250µl | 1-10µl                            | 5-120µl | 10-300µl | 50-1000µl | 0.2-10µl                         | 5-120µl | 10-300µl | 0.5-10µl                          | 2-20µl | 10-100µl | 20-200µl                        | 100-1000µl | 0.5-10µl | 5-100µl | 30-300µl |
| Volume        | Product No. | Universal Fit Filter Tips        |         |          |                                   |         |          |           |                                  |         |          |                                   |        |          |                                 |            |          |         |          |
| 0.5-10ul      | 20-0010     | ✓                                |         |          | ✓                                 |         |          |           | ✓                                |         |          |                                   |        |          |                                 |            | ✓        |         |          |
|               | 21-0010     | ✓                                |         |          | ✓                                 |         |          |           | ✓                                |         |          |                                   |        |          |                                 |            | ✓        |         |          |
| 1-200ul       | 20-0200     |                                  | ✓       | ✓        |                                   | ✓       | ✓        |           |                                  | ✓       | ✓        |                                   | ✓      | ✓        | ✓                               |            |          | ✓       | ✓        |
|               | 21-0200     |                                  | ✓       | ✓        |                                   | ✓       | ✓        |           |                                  | ✓       | ✓        |                                   | ✓      | ✓        | ✓                               |            |          | ✓       | ✓        |
| 100-1000ul    | 20-1000     |                                  |         |          |                                   |         |          | ✓         |                                  |         |          |                                   |        |          |                                 | ✓          |          |         |          |
|               | 21-1000     |                                  |         |          |                                   |         |          | ✓         |                                  |         |          |                                   |        |          |                                 | ✓          |          |         |          |
| 100-1250ul    | 20-1250     |                                  |         |          |                                   |         |          | ✓         |                                  |         |          |                                   |        |          |                                 | ✓          |          |         |          |
|               | 20-1251     |                                  |         |          |                                   |         |          | ✓         |                                  |         |          |                                   |        |          |                                 | ✓          |          |         |          |
|               | 21-1250     |                                  |         |          |                                   |         |          | ✓         |                                  |         |          |                                   |        |          |                                 | ✓          |          |         |          |

| Pipette Brand |             | Capp Aero (Multichannel)  |          |        |          |          |          | BrandTech Transferpette (Single channel) |         |          |        |        |          |          | MatrixImpact/Impact(s/mc) |            |          |         |          |
|---------------|-------------|---------------------------|----------|--------|----------|----------|----------|------------------------------------------|---------|----------|--------|--------|----------|----------|---------------------------|------------|----------|---------|----------|
| Model / Range |             | 0.2-2µl                   | 0.5-10µl | 5-50µl | 10-100µl | 25-200µl | 50-300µl | 100-1000µl                               | 0.1-1µl | 0.5-10µl | 2-20µl | 5-50µl | 10-100µl | 20-200µl | 25-250µl                  | 100-1000µl | 20-200µl | 5-100µl | 30-300µl |
| Volume        | Product No. | Universal Fit Filter Tips |          |        |          |          |          |                                          |         |          |        |        |          |          |                           |            |          |         |          |
| 0.5-10ul      | 20-0010     | ✓                         | ✓        |        |          |          |          |                                          | ✓       | ✓        |        |        |          |          |                           |            |          |         |          |
|               | 21-0010     | ✓                         | ✓        |        |          |          |          |                                          | ✓       | ✓        |        |        |          |          |                           |            |          |         |          |
| 1-200ul       | 20-0200     |                           |          | ✓      | ✓        | ✓        | ✓        |                                          |         |          |        | ✓      | ✓        | ✓        | ✓                         |            | ✓        |         |          |
|               | 21-0200     |                           |          | ✓      | ✓        | ✓        | ✓        |                                          |         |          |        | ✓      | ✓        | ✓        | ✓                         |            | ✓        |         |          |
| 100-1000ul    | 20-1000     |                           |          |        |          |          |          | ✓                                        |         |          |        |        |          |          |                           | ✓          |          |         |          |
|               | 21-1000     |                           |          |        |          |          |          | ✓                                        |         |          |        |        |          |          |                           | ✓          |          |         |          |
| 100-1250ul    | 20-1250     |                           |          |        |          |          |          | ✓                                        |         |          |        |        |          |          |                           | ✓          |          |         | ✓        |
|               | 20-1251     |                           |          |        |          |          |          | ✓                                        |         |          |        |        |          |          |                           | ✓          |          |         | ✓        |
|               | 21-1250     |                           |          |        |          |          |          | ✓                                        |         |          |        |        |          |          |                           | ✓          |          |         | ✓        |

# Pipet Tip Compatibility

| Pipette Brand |             | BrandTech Transferpette<br>(Multichannel) |        |          |        |          |          | Rainin EDP3<br>(Single channel) |                       |                    |                     |                     |                        | Rainin Latch Mode |                  |                |                   |                   |                      |
|---------------|-------------|-------------------------------------------|--------|----------|--------|----------|----------|---------------------------------|-----------------------|--------------------|---------------------|---------------------|------------------------|-------------------|------------------|----------------|-------------------|-------------------|----------------------|
| Model / Range |             | 0.5-10µl                                  | 2-20µl | 2.5-25µl | 5-50µl | 10-100µl | 20-200µl | 30-300µl                        | E-MIC-10,<br>0.5-10µl | E2-25,<br>2.5-25µl | E2-100,<br>10-100µl | E2-250,<br>25-250µl | E2-1000,<br>100-1000µl | R2,<br>0.1-2µl    | R10,<br>0.2-10µl | R20,<br>2-20µl | R100,<br>10-100µl | R200,<br>20-200µl | R1000,<br>100-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips                 |        |          |        |          |          |                                 |                       |                    |                     |                     |                        |                   |                  |                |                   |                   |                      |
| 0.5-10ul      | 20-0010     | ✓                                         | ✓      |          |        |          |          |                                 | ✓                     |                    |                     |                     |                        | ✓                 | ✓                |                |                   |                   |                      |
|               | 21-0010     | ✓                                         | ✓      |          |        |          |          |                                 | ✓                     |                    |                     |                     |                        | ✓                 | ✓                |                |                   |                   |                      |
| 1-200ul       | 20-0200     |                                           |        | ✓        | ✓      | ✓        | ✓        |                                 |                       | ✓                  | ✓                   | ✓                   |                        |                   |                  | ✓              | ✓                 | ✓                 |                      |
|               | 21-0200     |                                           |        | ✓        | ✓      | ✓        | ✓        |                                 |                       | ✓                  | ✓                   | ✓                   |                        |                   |                  | ✓              | ✓                 | ✓                 |                      |
| 100-1000ul    | 20-1000     |                                           |        |          |        |          |          | ✓                               |                       |                    |                     |                     | ✓                      |                   |                  |                |                   |                   | ✓                    |
|               | 21-1000     |                                           |        |          |        |          |          | ✓                               |                       |                    |                     |                     | ✓                      |                   |                  |                |                   |                   | ✓                    |
| 100-1250ul    | 20-1250     |                                           |        |          |        |          |          | ✓                               |                       |                    |                     |                     | ✓                      |                   |                  |                |                   |                   | ✓                    |
|               | 20-1251     |                                           |        |          |        |          |          | ✓                               |                       |                    |                     |                     | ✓                      |                   |                  |                |                   |                   | ✓                    |
|               | 21-1250     |                                           |        |          |        |          |          | ✓                               |                       |                    |                     |                     | ✓                      |                   |                  |                |                   |                   | ✓                    |

| Pipette Brand |             | Labnet Labpette (Single channel) |          |        |          |          |            | Labnet Labpette (Multichannel) |        |          |          | Labnet Biopette E (Single channel) |        |          |            | Labnet Biopette E (Multichannel) |        |          |            |
|---------------|-------------|----------------------------------|----------|--------|----------|----------|------------|--------------------------------|--------|----------|----------|------------------------------------|--------|----------|------------|----------------------------------|--------|----------|------------|
| Model / Range |             | 0.1-2µl                          | 0.5-10µl | 2-20µl | 10-100µl | 20-200µl | 100-1000µl | 1-10µl                         | 5-50µl | 20-200µl | 50-300µl | 0.5-10µl                           | 2-20µl | 10-200µl | 100-1000µl | 0.5-10µl                         | 2-20µl | 10-200µl | 100-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips        |          |        |          |          |            |                                |        |          |          |                                    |        |          |            |                                  |        |          |            |
| 0.5-10ul      | 20-0010     | ✓                                | ✓        |        |          |          |            | ✓                              |        |          |          | ✓                                  | ✓      |          |            | ✓                                | ✓      |          |            |
|               | 21-0010     | ✓                                | ✓        |        |          |          |            | ✓                              |        |          |          | ✓                                  | ✓      |          |            | ✓                                | ✓      |          |            |
| 1-200ul       | 20-0200     |                                  |          | ✓      | ✓        | ✓        |            |                                | ✓      | ✓        | ✓        |                                    |        | ✓        |            |                                  |        | ✓        |            |
|               | 21-0200     |                                  |          | ✓      | ✓        | ✓        |            |                                | ✓      | ✓        | ✓        |                                    |        | ✓        |            |                                  |        | ✓        |            |
| 100-1000ul    | 20-1000     |                                  |          |        |          |          | ✓          |                                |        |          |          |                                    |        |          | ✓          |                                  |        |          | ✓          |
|               | 21-1000     |                                  |          |        |          |          | ✓          |                                |        |          |          |                                    |        |          | ✓          |                                  |        |          | ✓          |
| 100-1250ul    | 20-1250     |                                  |          |        |          |          | ✓          |                                |        |          |          |                                    |        |          | ✓          |                                  |        |          | ✓          |
|               | 20-1251     |                                  |          |        |          |          | ✓          |                                |        |          |          |                                    |        |          | ✓          |                                  |        |          | ✓          |
|               | 21-1250     |                                  |          |        |          |          | ✓          |                                |        |          |          |                                    |        |          | ✓          |                                  |        |          | ✓          |

| Pipette Brand |             | Eppendorf Electronic      |         |         |           | Socorex Calibra/Accura (Single channel) |         |          |        |        |          |          | VWR AND VWR Signature(Multichannel) |        |        |          | Rainin EDP3 (Multichannel) |                    |                    |                   |
|---------------|-------------|---------------------------|---------|---------|-----------|-----------------------------------------|---------|----------|--------|--------|----------|----------|-------------------------------------|--------|--------|----------|----------------------------|--------------------|--------------------|-------------------|
| Model / Range |             | 0.5-10µl                  | 5-100µl | 2-200µl | 50-1000µl | 0.1-2µl                                 | 0.2-2µl | 0.5-10µl | 1-10µl | 2-20µl | 10-100µl | 20-200µl | 100-1000µl                          | 1-10µl | 5-50µl | 20-200µl | 50-300µl                   | EP-M8-10, 0.5-10µl | EP-M8250, 10-250µl | P200-M8, 20-200µl |
| Volume        | Product No. | Universal Fit Filter Tips |         |         |           |                                         |         |          |        |        |          |          |                                     |        |        |          |                            |                    |                    |                   |
| 0.5-10ul      | 20-0010     | ✓                         |         |         |           | ✓                                       | ✓       | ✓        | ✓      |        |          |          |                                     | ✓      |        |          |                            | ✓                  |                    |                   |
|               | 21-0010     | ✓                         |         |         |           | ✓                                       | ✓       | ✓        | ✓      |        |          |          |                                     | ✓      |        |          |                            | ✓                  |                    |                   |
| 1-200ul       | 20-0200     |                           | ✓       | ✓       |           |                                         |         |          |        | ✓      | ✓        | ✓        |                                     |        | ✓      | ✓        | ✓                          |                    | ✓                  | ✓                 |
|               | 21-0200     |                           | ✓       | ✓       |           |                                         |         |          |        | ✓      | ✓        | ✓        |                                     |        | ✓      | ✓        | ✓                          |                    | ✓                  | ✓                 |
| 100-1000ul    | 20-1000     |                           |         |         | ✓         |                                         |         |          |        |        |          |          | ✓                                   |        |        |          |                            |                    |                    | ✓                 |
|               | 21-1000     |                           |         |         | ✓         |                                         |         |          |        |        |          |          | ✓                                   |        |        |          |                            |                    |                    | ✓                 |
| 100-1250ul    | 20-1250     |                           |         |         | ✓         |                                         |         |          |        |        |          |          | ✓                                   |        |        |          |                            |                    |                    |                   |
|               | 20-1251     |                           |         |         | ✓         |                                         |         |          |        |        |          |          | ✓                                   |        |        |          |                            |                    |                    |                   |
|               | 21-1250     |                           |         |         | ✓         |                                         |         |          |        |        |          |          | ✓                                   |        |        |          |                            |                    |                    |                   |

# Pipet Tip Compatibility

| Pipette Brand |             | VWR AND VWR Signature<br>(Single channel) |          |        |        |          |          |          |            | Socorex Calibra/Accura<br>(Multichaaanel) |        |          |          |          | EppendorfReference<br>(single channel) (multichannel) |        |          |          |            |
|---------------|-------------|-------------------------------------------|----------|--------|--------|----------|----------|----------|------------|-------------------------------------------|--------|----------|----------|----------|-------------------------------------------------------|--------|----------|----------|------------|
| Model / Range |             | 0.1-2µl                                   | 0.5-10µl | 2-20µl | 5-50µl | 10-100µl | 20-200µl | 50-250µl | 100-1000µl | 1-10µl                                    | 5-50µl | 10-100µl | 20-200µl | 50-200µl | 0.1-10µl                                              | 2-20µl | 10-100µl | 20-200µl | 100-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips                 |          |        |        |          |          |          |            |                                           |        |          |          |          |                                                       |        |          |          |            |
| 0.5-10ul      | 20-0010     | ✓                                         | ✓        |        |        |          |          |          |            | ✓                                         |        |          |          |          | ✓                                                     |        |          |          |            |
|               | 21-0010     | ✓                                         | ✓        |        |        |          |          |          |            | ✓                                         |        |          |          |          | ✓                                                     |        |          |          |            |
| 1-200ul       | 20-0200     |                                           |          | ✓      | ✓      | ✓        | ✓        | ✓        |            |                                           | ✓      | ✓        | ✓        | ✓        |                                                       | ✓      | ✓        | ✓        |            |
|               | 21-0200     |                                           |          | ✓      | ✓      | ✓        | ✓        | ✓        |            |                                           | ✓      | ✓        | ✓        | ✓        |                                                       | ✓      | ✓        | ✓        |            |
| 100-1000ul    | 20-1000     |                                           |          |        |        |          |          |          | ✓          |                                           |        |          |          |          |                                                       | ✓      | ✓        | ✓        |            |
|               | 21-1000     |                                           |          |        |        |          |          |          | ✓          |                                           |        |          |          |          |                                                       |        |          |          | ✓          |
| 100-1250ul    | 20-1250     |                                           |          |        |        |          |          |          | ✓          |                                           |        |          |          |          |                                                       |        |          |          | ✓          |
|               | 20-1251     |                                           |          |        |        |          |          |          | ✓          |                                           |        |          |          |          |                                                       |        |          |          | ✓          |
|               | 21-1250     |                                           |          |        |        |          |          |          | ✓          |                                           |        |          |          |          |                                                       |        |          |          | ✓          |

| Pipette Brand |             | Eppendorf<br>Research (singlechannel) |          |        |          |          |            | Eppendorf<br>Research (multichannel) |         |          |          |          |           | Thermo Scientif Finnpiptette<br>(single channel) |          |        |          |            |            |
|---------------|-------------|---------------------------------------|----------|--------|----------|----------|------------|--------------------------------------|---------|----------|----------|----------|-----------|--------------------------------------------------|----------|--------|----------|------------|------------|
| Model / Range |             | 0.1-2.5µl                             | 0.5-10µl | 2-20µl | 10-100µl | 25-200µl | 100-1000µl | 0.5-10µl                             | 5-100µl | 50-100µl | 20-300µl | 30-300µl | 50-1200µl | 0.2-10µl                                         | 0.5-10µl | 5-40µl | 40-200µl | 100-1000µl | 200-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips             |          |        |          |          |            |                                      |         |          |          |          |           |                                                  |          |        |          |            |            |
| 0.5-10ul      | 20-0010     | ✓                                     | ✓        |        |          |          |            | ✓                                    |         |          |          |          |           | ✓                                                | ✓        |        |          |            |            |
|               | 21-0010     | ✓                                     | ✓        |        |          |          |            | ✓                                    |         |          |          |          |           | ✓                                                | ✓        |        |          |            |            |
| 1-200ul       | 20-0200     |                                       |          |        | ✓        | ✓        |            |                                      |         | ✓        | ✓        | ✓        | ✓         |                                                  |          | ✓      | ✓        |            |            |
|               | 21-0200     |                                       |          |        | ✓        | ✓        |            |                                      |         | ✓        | ✓        | ✓        | ✓         |                                                  |          | ✓      | ✓        |            |            |
| 100-1000ul    | 20-1000     |                                       |          |        | ✓        | ✓        |            |                                      |         | ✓        | ✓        | ✓        | ✓         |                                                  |          | ✓      | ✓        |            |            |
|               | 21-1000     |                                       |          |        |          |          | ✓          |                                      |         |          |          |          | ✓         |                                                  |          |        |          | ✓          | ✓          |
| 100-1250ul    | 20-1250     |                                       |          |        |          |          | ✓          |                                      |         |          |          |          | ✓         |                                                  |          |        |          | ✓          | ✓          |
|               | 20-1251     |                                       |          |        |          |          | ✓          |                                      |         |          |          |          |           |                                                  |          |        |          | ✓          | ✓          |
|               | 21-1250     |                                       |          |        |          |          | ✓          |                                      |         |          |          |          |           |                                                  |          |        |          | ✓          | ✓          |

| Pipette Brand |             | Thermo Scientif<br>Finnpiptette (multichannel) |          | Gilson Pipetma (single channel) |               |             |                |                |                   |
|---------------|-------------|------------------------------------------------|----------|---------------------------------|---------------|-------------|----------------|----------------|-------------------|
| Model / Range |             | 5-50µl                                         | 50-300µl | P20, 0.1-2µl                    | P10, 0.2-10µl | P20, 2-20µl | P100, 10-100µl | P200, 20-200µl | P1000, 100-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips                      |          |                                 |               |             |                |                |                   |
| 0.5-10ul      | 20-0010     |                                                |          | ✓                               | ✓             |             |                |                |                   |
|               | 21-0010     |                                                |          | ✓                               | ✓             |             |                |                |                   |
| 1-200ul       | 20-0200     | ✓                                              | ✓        |                                 |               | ✓           | ✓              | ✓              |                   |
|               | 21-0200     | ✓                                              | ✓        |                                 |               | ✓           | ✓              | ✓              |                   |
| 100-1000ul    | 20-1000     | ✓                                              | ✓        |                                 |               | ✓           | ✓              | ✓              |                   |
|               | 21-1000     |                                                |          |                                 |               |             |                |                | ✓                 |
| 100-1250ul    | 20-1250     |                                                |          |                                 |               |             |                |                | ✓                 |
|               | 20-1251     |                                                |          |                                 |               |             |                |                | ✓                 |
|               | 21-1250     |                                                |          |                                 |               |             |                |                | ✓                 |

# Pipet Tip Compatibility

| Pipette Brand |             | Eppendorf Research Plus   |          |        |          |          |            | Gilson Pipetman(PG) |             |             |                |                |                   | Gilson Pipetman(PN) |             |             |                |                |                   | Gilson Pipetman(PL) |             |             |                |                |                   |
|---------------|-------------|---------------------------|----------|--------|----------|----------|------------|---------------------|-------------|-------------|----------------|----------------|-------------------|---------------------|-------------|-------------|----------------|----------------|-------------------|---------------------|-------------|-------------|----------------|----------------|-------------------|
| Model / Range |             | 0.1-2.5µl                 | 0.5-10µl | 2-20µl | 10-100µl | 20-200µl | 100-1000µl | P2, 0.2-2µl         | P10, 1-10µl | P20, 2-20µl | P100, 10-100µl | P200, 20-200µl | P1000, 100-1000µl | P2, 0.2-2µl         | P10, 1-10µl | P20, 2-20µl | P100, 10-100µl | P200, 20-200µl | P1000, 100-1000µl | P2, 0.2-2µl         | P10, 1-10µl | P20, 2-20µl | P100, 10-100µl | P200, 20-200µl | P1000, 100-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips |          |        |          |          |            |                     |             |             |                |                |                   |                     |             |             |                |                |                   |                     |             |             |                |                |                   |
| 0.5-10ul      | 22-0010S    | ✓                         | ✓        |        |          |          |            | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   |
|               | 22-0011S    | ✓                         | ✓        |        |          |          |            | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   |
| 1-20ul        | 22-0020S    |                           |          | ✓      |          |          |            |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |
|               | 22-0021S    |                           |          | ✓      |          |          |            |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |
| 1-50ul        | 22-0050S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
|               | 22-0051S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
| 1-100ul       | 22-0100S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
|               | 22-0101S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
| 1-300ul       | 22-0200S    |                           |          |        |          | ✓        |            |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |
|               | 22-0201S    |                           |          |        |          | ✓        |            |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |
| 100-1000ul    | 22-1000S    |                           |          |        |          |          | ✓          |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |
|               | 22-1001S    |                           |          |        |          |          | ✓          |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |

| Pipette Brand |             | Thermo Finnpiptette (F1)  |        |              |        |              |        |              |          |          |          | Thermo Finnpiptette (F2) |               |        |              |        |              |        |              |          |          | Thermo Finnpiptette (F3) |         |        |              |        |              |        |              |          |          |            |
|---------------|-------------|---------------------------|--------|--------------|--------|--------------|--------|--------------|----------|----------|----------|--------------------------|---------------|--------|--------------|--------|--------------|--------|--------------|----------|----------|--------------------------|---------|--------|--------------|--------|--------------|--------|--------------|----------|----------|------------|
| Model / Range |             | 0.2-2µl                   | 1-10µl | 1-10µl(µTip) | 2-20µl | 2-20µl(µTip) | 5-50µl | 5-50µl(µTip) | 10-100µl | 20-200µl | 30-300µl | 100-1000µl               | 0.2-2µl(µTip) | 1-10µl | 1-10µl(µTip) | 2-20µl | 2-20µl(µTip) | 5-50µl | 5-50µl(µTip) | 10-100µl | 20-200µl | 100-1000µl               | 0.2-2µl | 1-10µl | 1-10µl(µTip) | 2-20µl | 2-20µl(µTip) | 5-50µl | 5-50µl(µTip) | 10-100µl | 20-200µl | 100-1000µl |
| Volume        | Product No. | Universal Fit Filter Tips |        |              |        |              |        |              |          |          |          |                          |               |        |              |        |              |        |              |          |          |                          |         |        |              |        |              |        |              |          |          |            |
| 0.5-10ul      | 22-0010S    | ✓                         | ✓      |              |        |              |        |              |          |          |          |                          | ✓             |        |              |        |              |        |              |          |          |                          | ✓       | ✓      |              |        |              |        |              |          |          |            |
|               | 22-0011S    | ✓                         | ✓      |              |        |              |        |              |          |          |          |                          | ✓             |        |              |        |              |        |              |          |          |                          | ✓       | ✓      |              |        |              |        |              |          |          |            |
| 1-20ul        | 22-0020S    |                           |        |              | ✓      |              |        |              |          |          |          |                          |               |        |              | ✓      |              |        |              |          |          |                          |         |        |              | ✓      |              |        |              |          |          |            |
|               | 22-0021S    |                           |        |              | ✓      |              |        |              |          |          |          |                          |               |        |              | ✓      |              |        |              |          |          |                          |         |        | ✓            |        |              |        |              |          |          |            |
| 1-50ul        | 22-0050S    |                           |        |              |        | ✓            |        |              | ✓        | ✓        |          |                          |               |        |              |        | ✓            |        |              | ✓        | ✓        |                          |         |        |              |        | ✓            |        |              | ✓        | ✓        |            |
|               | 22-0051S    |                           |        |              |        | ✓            |        |              | ✓        | ✓        |          |                          |               |        |              |        | ✓            | ✓      |              | ✓        | ✓        |                          |         |        |              |        | ✓            | ✓      |              | ✓        | ✓        |            |
| 1-100ul       | 22-0100S    |                           |        |              |        | ✓            |        |              | ✓        | ✓        |          |                          |               |        |              |        | ✓            |        |              | ✓        | ✓        |                          |         |        |              |        | ✓            |        |              | ✓        | ✓        |            |
|               | 22-0101S    |                           |        |              |        | ✓            |        |              | ✓        | ✓        |          |                          |               |        |              |        | ✓            |        |              | ✓        | ✓        |                          |         |        |              |        | ✓            |        |              | ✓        | ✓        |            |
| 1-300ul       | 22-0200S    |                           |        |              |        |              |        |              |          |          | ✓        |                          |               |        |              |        |              |        |              |          | ✓        | ✓                        |         |        |              |        |              |        |              | ✓        | ✓        |            |
|               | 22-0201S    |                           |        |              |        |              |        |              |          |          | ✓        |                          |               |        |              |        |              |        |              |          | ✓        |                          |         |        |              |        |              |        |              | ✓        | ✓        |            |
| 100-1000ul    | 22-1000S    |                           |        |              |        |              |        |              |          |          |          | ✓                        |               |        |              |        |              |        |              |          |          | ✓                        |         |        |              |        |              |        |              |          |          | ✓          |
|               | 22-1001S    |                           |        |              |        |              |        |              |          |          |          | ✓                        |               |        |              |        |              |        |              |          |          | ✓                        |         |        |              |        |              |        |              |          |          | ✓          |

# Pipet Tip Compatibility

| Pipette Brand |             | Thermo Finnpiipette (F1)  |          |        |        |          |          |          |            | Thermo Finnpiipette (F2) |          |        |        |          |          |          |            | Thermo Finnpiipette (F3) |           |          |        |        |          |          |          |            |   |
|---------------|-------------|---------------------------|----------|--------|--------|----------|----------|----------|------------|--------------------------|----------|--------|--------|----------|----------|----------|------------|--------------------------|-----------|----------|--------|--------|----------|----------|----------|------------|---|
| Model / Range |             | 0.1-2.5µl                 | 0.5-10µl | 2-20µl | 5-50µl | 10-100µl | 20-200µl | 50-300µl | 100-1000µl | 0.1-3µl                  | 0.5-10µl | 2-20µl | 5-50µl | 10-100µl | 20-200µl | 30-300µl | 100-1000µl | 0.1-1µl                  | 0.1-2.5µl | 0.5-10µl | 2-20µl | 5-50µl | 10-100µl | 20-200µl | 30-300µl | 100-1000µl |   |
| Volume        | Product No. | Universal Fit Filter Tips |          |        |        |          |          |          |            |                          |          |        |        |          |          |          |            |                          |           |          |        |        |          |          |          |            |   |
| 0.5-10ul      | 22-0010S    | ✓                         | ✓        |        |        |          |          |          |            | ✓                        | ✓        |        |        |          |          |          |            |                          | ✓         | ✓        | ✓      |        |          |          |          |            |   |
|               | 22-0011S    | ✓                         | ✓        |        |        |          |          |          |            | ✓                        | ✓        |        |        |          |          |          |            |                          | ✓         | ✓        | ✓      |        |          |          |          |            |   |
| 1-20ul        | 22-0020S    |                           |          | ✓      |        |          |          |          |            |                          |          | ✓      |        |          |          |          |            |                          |           | ✓        |        |        |          |          |          |            |   |
|               | 22-0021S    |                           |          | ✓      |        |          |          |          |            |                          |          | ✓      |        |          |          |          |            |                          |           | ✓        |        |        |          |          |          |            |   |
| 1-50ul        | 22-0050S    |                           |          |        | ✓      | ✓        | ✓        | ✓        |            |                          |          |        | ✓      | ✓        | ✓        | ✓        |            |                          |           |          |        | ✓      | ✓        | ✓        | ✓        |            |   |
|               | 22-0051S    |                           |          |        | ✓      | ✓        | ✓        | ✓        |            |                          |          |        | ✓      | ✓        | ✓        | ✓        |            |                          |           |          |        | ✓      | ✓        | ✓        | ✓        |            |   |
| 1-100ul       | 22-0100S    |                           |          |        |        | ✓        | ✓        | ✓        |            |                          |          |        |        | ✓        | ✓        | ✓        |            |                          |           |          |        |        | ✓        | ✓        | ✓        |            |   |
|               | 22-0101S    |                           |          |        |        | ✓        | ✓        | ✓        |            |                          |          |        |        | ✓        | ✓        | ✓        |            |                          |           |          |        |        | ✓        | ✓        | ✓        |            |   |
| 1-300ul       | 22-0200S    |                           |          |        |        |          | ✓        | ✓        |            |                          |          |        |        |          | ✓        | ✓        |            |                          |           |          |        |        |          | ✓        | ✓        |            |   |
|               | 22-0201S    |                           |          |        |        |          | ✓        | ✓        |            |                          |          |        |        |          | ✓        | ✓        |            |                          |           |          |        |        |          | ✓        | ✓        |            |   |
| 100-1000ul    | 22-1000S    |                           |          |        |        |          |          |          | ✓          |                          |          |        |        |          |          |          | ✓          |                          |           |          |        |        |          |          |          |            | ✓ |
|               | 22-1001S    |                           |          |        |        |          |          |          | ✓          |                          |          |        |        |          |          |          | ✓          |                          |           |          |        |        |          |          |          | ✓          |   |

| Pipette Brand |             | Eppendorf Research Plus   |          |        |          |          |            | Gilson Pipetman(PG) |             |             |                |                |                   | Gilson Pipetman(PN) |             |             |                |                |                   | Gilson Pipetman(PL) |             |             |                |                |                   |
|---------------|-------------|---------------------------|----------|--------|----------|----------|------------|---------------------|-------------|-------------|----------------|----------------|-------------------|---------------------|-------------|-------------|----------------|----------------|-------------------|---------------------|-------------|-------------|----------------|----------------|-------------------|
| Model / Range |             | 0.1-2.5ul                 | 0.5-10ul | 2-20ul | 10-100ul | 20-200ul | 100-1000ul | P2, 0.2-2ul         | P10, 1-10ul | P20, 2-20ul | P100, 10-100ul | P200, 20-200ul | P1000, 100-1000ul | P2, 0.2-2ul         | P10, 1-10ul | P20, 2-20ul | P100, 10-100ul | P200, 20-200ul | P1000, 100-1000ul | P2, 0.2-2ul         | P10, 1-10ul | P20, 2-20ul | P100, 10-100ul | P200, 20-200ul | P1000, 100-1000ul |
| Volume        | Product No. | Universal Fit Filter Tips |          |        |          |          |            |                     |             |             |                |                |                   |                     |             |             |                |                |                   |                     |             |             |                |                |                   |
| 0.5-10ul      | 23-0010S    | ✓                         | ✓        |        |          |          |            | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   |
|               | 23-0011S    | ✓                         | ✓        |        |          |          |            | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   | ✓                   | ✓           |             |                |                |                   |
| 1-20ul        | 23-0020S    |                           |          | ✓      |          |          |            |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |
|               | 23-0021S    |                           |          | ✓      |          |          |            |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |                     |             | ✓           |                |                |                   |
| 1-50ul        | 23-0050S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
|               | 23-0051S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
| 1-100ul       | 23-0100S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
|               | 23-0101S    |                           |          |        | ✓        |          |            |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |                     |             |             | ✓              | ✓              |                   |
| 1-300ul       | 23-0200S    |                           |          |        |          | ✓        |            |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |
|               | 23-0201S    |                           |          |        |          | ✓        |            |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |                     |             |             |                | ✓              |                   |
| 100-1000ul    | 23-1000S    |                           |          |        |          |          | ✓          |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |
|               | 23-1001S    |                           |          |        |          |          | ✓          |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |                     |             |             |                |                | ✓                 |

# Pipet Tip Compatibility

| Pipette Brand |             | Thermo Finnpiquette (F1)  |        |               |        |               |        |               |          |          |          | Thermo Finnpiquette (F2) |                |        |               |        |               |        |               |          |          | Thermo Finnpiquette (F3) |         |        |               |        |               |        |               |          |          |            |
|---------------|-------------|---------------------------|--------|---------------|--------|---------------|--------|---------------|----------|----------|----------|--------------------------|----------------|--------|---------------|--------|---------------|--------|---------------|----------|----------|--------------------------|---------|--------|---------------|--------|---------------|--------|---------------|----------|----------|------------|
| Model / Range |             | 0.2-2ul                   | 1-10ul | 1-10ul(μTip ) | 2-20ul | 2-20ul(μTip ) | 5-50ul | 5-50ul(μTip ) | 10-100ul | 20-200ul | 30-300ul | 100-1000ul               | 0.2-2ul(μTip ) | 1-10ul | 1-10ul(μTip ) | 2-20ul | 2-20ul(μTip ) | 5-50ul | 5-50ul(μTip ) | 10-100ul | 20-200ul | 100-1000ul               | 0.2-2ul | 1-10ul | 1-10ul(μTip ) | 2-20ul | 2-20ul(μTip ) | 5-50ul | 5-50ul(μTip ) | 10-100ul | 20-200ul | 100-1000ul |
|               |             | Universal Fit Filter Tips |        |               |        |               |        |               |          |          |          |                          |                |        |               |        |               |        |               |          |          |                          |         |        |               |        |               |        |               |          |          |            |
| Volume        | Product No. |                           |        |               |        |               |        |               |          |          |          |                          |                |        |               |        |               |        |               |          |          |                          |         |        |               |        |               |        |               |          |          |            |
|               | 23-0010S    | ✓                         | ✓      |               |        |               |        |               |          |          |          |                          |                | ✓      |               |        |               |        |               |          |          |                          |         | ✓      | ✓             |        |               |        |               |          |          |            |
| 0.5-10ul      | 23-0011S    | ✓                         | ✓      |               |        |               |        |               |          |          |          |                          | ✓              |        |               |        |               |        |               |          |          |                          |         | ✓      | ✓             |        |               |        |               |          |          |            |
|               | 23-0020S    |                           |        |               |        | ✓             |        |               |          |          |          |                          |                |        |               |        |               |        |               |          |          |                          |         |        |               |        | ✓             |        |               |          |          |            |
| 1-20ul        | 23-0021S    |                           |        |               | ✓      |               |        |               |          |          |          |                          |                |        |               | ✓      |               |        |               |          |          |                          |         |        |               | ✓      |               |        |               |          |          |            |
|               | 23-0050S    |                           |        |               |        |               | ✓      |               | ✓        | ✓        |          |                          |                |        |               |        |               | ✓      |               |          |          |                          |         |        |               |        |               | ✓      |               | ✓        | ✓        |            |
| 1-50ul        | 23-0051S    |                           |        |               |        | ✓             |        | ✓             | ✓        | ✓        |          |                          |                |        |               |        | ✓             |        |               |          |          |                          |         |        |               |        |               | ✓      |               | ✓        | ✓        |            |
|               | 23-0100S    |                           |        |               |        | ✓             |        | ✓             | ✓        | ✓        |          |                          |                |        |               |        | ✓             |        |               |          |          |                          |         |        |               |        |               | ✓      |               | ✓        | ✓        |            |
| 1-100ul       | 23-0101S    |                           |        |               |        | ✓             |        | ✓             | ✓        | ✓        |          |                          |                |        |               |        | ✓             |        |               |          |          |                          |         |        |               |        |               | ✓      |               | ✓        | ✓        |            |
|               | 23-0200S    |                           |        |               |        |               |        |               |          |          | ✓        |                          |                |        |               |        |               |        |               |          |          |                          |         |        |               |        |               |        |               | ✓        | ✓        |            |
| 1-300ul       | 23-0201S    |                           |        |               |        |               |        |               |          |          | ✓        |                          |                |        |               |        |               |        |               |          |          |                          |         |        |               |        |               |        |               | ✓        | ✓        |            |
|               | 23-1000S    |                           |        |               |        |               |        |               |          |          |          | ✓                        |                |        |               |        |               |        |               |          |          |                          |         |        |               |        |               |        |               |          | ✓        |            |
| 100-1000ul    | 23-1001S    |                           |        |               |        |               |        |               |          |          |          | ✓                        |                |        |               |        |               |        |               |          |          |                          |         |        |               |        |               |        |               |          | ✓        |            |

| Pipette Brand |             | BIOHIT Proline            |          |        |        |          |          |          |            | BIOHIT Proline (Plus) |          |        |        |          |          |          |            | BRAND   |           |          |        |        |          |          |          |            |  |
|---------------|-------------|---------------------------|----------|--------|--------|----------|----------|----------|------------|-----------------------|----------|--------|--------|----------|----------|----------|------------|---------|-----------|----------|--------|--------|----------|----------|----------|------------|--|
| Model / Range |             | 0.1-2.5ul                 | 0.5-10ul | 2-20ul | 5-50ul | 10-100ul | 20-200ul | 50-300ul | 100-1000ul | 0.1-3ul               | 0.5-10ul | 2-20ul | 5-50ul | 10-100ul | 20-200ul | 30-300ul | 100-1000ul | 0.1-1ul | 0.1-2.5ul | 0.5-10ul | 2-20ul | 5-50ul | 10-100ul | 20-200ul | 30-300ul | 100-1000ul |  |
| Volume        | Product No. | Universal Fit Filter Tips |          |        |        |          |          |          |            |                       |          |        |        |          |          |          |            |         |           |          |        |        |          |          |          |            |  |
| 0.5-10ul      | 23-0010S    | ✓                         | ✓        |        |        |          |          |          |            | ✓                     | ✓        |        |        |          |          |          |            | ✓       | ✓         | ✓        |        |        |          |          |          |            |  |
|               | 23-0011S    | ✓                         | ✓        |        |        |          |          |          |            | ✓                     | ✓        |        |        |          |          |          |            | ✓       | ✓         | ✓        |        |        |          |          |          |            |  |
| 1-20ul        | 23-0020S    |                           |          | ✓      |        |          |          |          |            |                       |          | ✓      |        |          |          |          |            |         |           |          | ✓      |        |          |          |          |            |  |
|               | 23-0021S    |                           |          | ✓      |        |          |          |          |            |                       |          | ✓      |        |          |          |          |            |         |           |          | ✓      |        |          |          |          |            |  |
| 1-50ul        | 23-0050S    |                           |          |        | ✓      | ✓        | ✓        | ✓        |            |                       |          |        | ✓      | ✓        | ✓        | ✓        |            |         |           |          |        | ✓      | ✓        | ✓        | ✓        |            |  |
|               | 23-0051S    |                           |          |        | ✓      | ✓        | ✓        | ✓        |            |                       |          |        | ✓      | ✓        | ✓        | ✓        |            |         |           |          |        | ✓      | ✓        | ✓        | ✓        |            |  |
| 1-100ul       | 23-0100S    |                           |          |        |        | ✓        | ✓        | ✓        |            |                       |          |        |        | ✓        | ✓        | ✓        |            |         |           |          |        |        | ✓        | ✓        | ✓        |            |  |
|               | 23-0101S    |                           |          |        |        | ✓        | ✓        | ✓        |            |                       |          |        |        | ✓        | ✓        | ✓        |            |         |           |          |        |        | ✓        | ✓        | ✓        |            |  |
| 1-300ul       | 23-0200S    |                           |          |        |        |          | ✓        | ✓        |            |                       |          |        |        |          | ✓        | ✓        |            |         |           |          |        |        |          | ✓        | ✓        |            |  |
|               | 23-0201S    |                           |          |        |        |          | ✓        | ✓        |            |                       |          |        |        |          | ✓        | ✓        |            |         |           |          |        |        |          | ✓        | ✓        |            |  |
| 100-1000ul    | 23-1000S    |                           |          |        |        |          |          |          | ✓          |                       |          |        |        |          |          |          | ✓          |         |           |          |        |        |          |          |          | ✓          |  |
|               | 23-1001S    |                           |          |        |        |          |          |          | ✓          |                       |          |        |        |          |          |          | ✓          |         |           |          |        |        |          |          |          | ✓          |  |