

Clean Trust™ Pipette Swab


The Pipette Swab is an innovative laboratory consumable that integrates both sampling and pipetting functions, specifically designed for pathogen sample collection. It is particularly suitable for environmental monitoring and microbial sampling in industries such as food and beverage.



Benefits

- **Precision** - Built-in graduations on the pipette allow for accurate and convenient liquid handling.
- **Efficiency** - Preloaded with diluent to simplify workflow and enhance sampling efficiency.
- **Sterility** - Irradiation-sterilized, ready to use straight out of the package for guaranteed aseptic sampling.
- **Viability** - Contains 10 mL of sterile buffer solution to repair bacteria and maintain sample activity.
- **Portability** - Easy to use, reduces manual steps, and lowers overall testing costs.
- **Versatility** - Compatible with petri dishes, test films, and more for rapid quantitative and qualitative microbial analysis.

Procedure



2 ~ 25°C, valid for 12 months.

Product No.	Description	Specification
CT-DE2510P	D/E Neutralizing Broth, 10mL	10 pcs/box
CT-BPW2510P	Buffered Peptone Water, 10mL	10 pcs/box
CT-NB2510P	Neutralizing Buffer, 10mL	10 pcs/box
CT-PBS2510P	Phosphate Buffered Saline, 10mL	10 pcs/box

LABTEK Services Ltd.

Tel: 0044 1935 473824

Email: sales@labtekservices.co.uk

108 Preston Grove, Yeovil, Somerset BA202DA

